


Achieve easy, fast and reliable  
particle/cell size analysis with  
NanoCuvette™ S

Introduction

# Before we begin



We will send you the  
recording



Submit your questions  
anytime. We'll do Q&A in  
the end



Contact us:  
[webinar@avantorsciences.com](mailto:webinar@avantorsciences.com)

## The speakers



Christopher Lüscher – Interim CEO, PhD



Sangita Khatri - Research Scientist




# Achieve easy, fast and reliable particle size analysis with NanoCuvette™ S

Upgrade your existing spectrophotometer with NanoCuvette™ S and cloud computing - our take on affordable particle and cell size analysis.



**cphnano.com**  
Hørmarken 2, 3520 Farum  
Denmark.  
Tel: +45 36 99 27 46  
info@cphnano.com



## OUTLINE

- About the company
- Absorption spectroscopy (Beer-Lambert law)
- VWR Collection spectrophotometers
- Particle Size Analysis
- NanoCuvette™ S
- SpectroWorks™
- Conclusion



**Sangita Khatri**

Research Scientist  
sakha@cphnano.com




**Christopher Lüscher, PhD**

Interim CEO  
chrлу@cphnano.com



## EASY, FAST AND RELIABLE PARTICLE SIZE ANALYSIS




Poll Question 1

**Do you have a cuvette-based UV-Vis spectrophotometer in your lab?**

## ABOUT THE COMPANY




- **Founded in 2015, cphnano is a Danish labtech company** from Copenhagen that develops digital laboratory analysis and next-generation UV-Vis for the smart lab of the future.
- **We work with upgrading UV-Vis spectrophotometry** across all brands to be at the forefront at the technology development for the smart lab of the future and Industry 4.0.
- **We sell consumables and digital upgrades** for laboratory equipment through our trusted resellers and partners for an integrated, seamless digital experience.





## TRUSTED BY



Customers and partners of cphnano range from biotech and life science companies to manufacturers, R&D labs and SME's as well as education institutions and university laboratories.

## PRODUCTS AND SERVICES



### NanoCuvette™ One and S (10 pcs.)

- Commodities
- Unique ID
- 8.5 mm and 15/20 mm versions
- Including SpectroWorks™ full functionality

### SpectroWorks™

- Online software for spectrophotometry
- License included with NanoCuvette™ purchase
- Supports basic plastic, quartz or glass cuvettes as well
- Includes instrument light simulations for better performance
- Drag and drop easy
- All future updates included

### SpectroLink™

- Instrument control from smartphone, tablet or computer.
- Best possible data quality.
- No installation.



## NanoCuvette™ S




## SpectroWorks™



## THE SPECTROPHOTOMETER BUILT FOR MEASURING BEER-LAMBERT LAW

- Spectrophotometry is one of the most widely used analytical procedures in many labs.
- It is used to estimate the level of a substance in a solution and is ideal for simple routine determination of small quantities of materials such as DNA, RNA, proteins, enzymes and many other substances.





## ABSORPTION SPECTROSCOPY AND THE BEER-LAMBERT LAW

The Beer-Lambert law relates the attenuation of light to the properties of the material through which the light is travelling:

$$A = \epsilon c l + A_0$$

where  $A$  is the absorbance,  $\epsilon$  is molar extinction coefficient,  $c$  is the concentration,  $l$  is the path length and  $A_0$  is the background absorbance.



# VWR SPECTROPHOTOMETER COLLECTION LINE-UP



| Model                            | V-1200   | UV1600-PC | 3100-PC  | 6300-PC  | M4       | P4       | PV4      |
|----------------------------------|----------|-----------|----------|----------|----------|----------|----------|
| Type                             | Scanning | Scanning  | Scanning | Scanning | Scanning | Scanning | Scanning |
| Wavelengths (nm)                 | 325-1000 | 190-1100  | 190-1100 | 190-1100 | 190-1100 | 190-1100 | 190-1100 |
| Beam height (mm)                 | 15       | 15        | 15       | 15       | 15       | 15       | 15       |
| Minimum spectral band width (nm) | 4        | 4         | 2.0      | 1        | 2.0      | 2.0      | 2.0      |
| Wavelength accuracy (nm)         | 2.0      | 0.5       | 0.5      | 0.3      | 0.5      | 0.5      | 0.5      |
| Max scanning speed (nm/min)      | 4200     | 4200      | 3000     | 3000     |          | 4200     | 4200     |
| Connectivity interface           | USB      | USB       | USB      | USB      | USB      | USB      | USB      |
| NanoCuvette™ compatible          | ✓        | ✓         | ✓        | ✓        | ✓        | ✓        | ✓        |

## V-1200 / UV-1600PC

The V-1200 and UV-1600PC are basic Visible and UV/Vis spectrophotometers for the routine daily quality control requirements, in process control and teaching laboratories

## 3100-PC

The 3100-PC for advanced measurements in quality control, biochemical and pharmaceutical applications.

## 6300-PC

The 6300-PC with double-beam optical system, designed for pharmaceutical, biochemical and clinical use, new material development and quality control.

## M4

The M4 single split-beam spectrophotometer, with flashing xenon lamp, is an easy to use, highly reliable and accurate product.

## P4

The P4 UV/Visible spectrophotometer is an easy to use instrument with advanced performance

## PV4


The PV4 visible spectrophotometer is an easy to use instrument with advanced performance.

Poll Question 2


**Do you have particle size analysis  
equipment in your lab?**

## INTRODUCTION TO PARTICLE SIZE ANALYSIS

- Particle Size Analysis is the technical procedure to determine the list of values or a mathematical function that defines the relative amounts of particles dispersed in fluid sorted according to the size.
- Particle size influence the dissolution, absorption, stability, physical, chemical and pharmacokinetic properties of particles.
- Particle Size Analyzers (PSA) are based on different technologies including analysis of Brownian motion, high-definition image processing and Rayleigh and Mie scattering.



## PARTICLE SIZE ANALYSIS METHODS





Poll Question 3


**Can your instrument measure both  
the size and concentration at the  
same time?**

One solution for particle size and concentration analysis

**NanoCuvette™ S**


## THE NANOCUVETTE UNLOCKS THE FULL POTENTIAL OF THE SPECTROPHOTOMETER

**Spectrophotometry** is one of the most widely used analytical procedures in biochemistry.




**Adding an optical filter to the cuvette** upgrades the analytical capability of the instrument:

- **More complementary information on same sample.**
- **Easy and fast measurement of particle size and concentration.** Direct label-free size and concentration detection from 10 nm -10  $\mu\text{m}$ .



# MEASUREMENT PRINCIPLE – PHOTONIC CRYSTALS AND CLOUD COMPUTING

FIGURE 4: LIGHT SCATTERING





Particle and cell sizes can be measured, since light scattering depends on angle and wavelength.

**Dynamic Light Scattering (DLS)** measures scattering over time at given wavelength.

**Multi Angle Dynamic Light Scattering (MA-DLS)** measures scattering over time at given wavelength for different angles (Malvern ZetaSizer Ultra).

**NanoCuvette™ S** measures scattering at different wavelengths in different angles at given time.


FIGURE 6: MEASURED MILK SPECTRUM



- +500 million light simulations in SpectroWorks™.
- Rayleigh scattering (nanoparticles) and Mie scattering (microparticles) simulations accessed quickly with modern cloud computing architecture.
- Side A, side B and side D SLS of the cuvette.
- +500 combinations of particle and fluid materials.

## NANOCUVETTE™ S BENEFITS


- **Less sample needed for analysis.** 50  $\mu$ L - 200  $\mu$ L sample volume.
- **Easy and fast measurement of size and concentration.** Direct label-free size and concentration detection from 10 nm -10  $\mu$ m (Limit of detection is below 0.003 % particle concentration determination).
- **Replacing expensive equipment with just the spectrophotometer.**
- **Includes all basic cuvette features.**



Poll Question 4

**Which sample do you usually work with for particle size analysis?**

# VALIDATION OF NANOCUVETTE™ S BY PARTICLE ANALYTICAL



EXPERTS IN SIZE AND CRYSTALS

Copenhagen Nanosystems ApS

## TESTING AND VALIDATION OF NANOCUVETTE™ S, SPECTROWORKS™ TOGETHER WITH UV-VIS SPECTROPHOTOMETER

30-Jul-2021

PARTICLE ANALYTICAL ApS | Agerm Allé 3 | DK-2970 Hørsholm | VAT No.: DK25194934  
T: +45 4576 3060 | [www.particle.dk](http://www.particle.dk) | [info@particle.dk](mailto:info@particle.dk)

Page 1 of 9



**PARTICLE ANALYTICAL**  
EXPERTS IN SIZE AND CRYSTALS

TESTING AND VALIDATION OF NANOCUVETTE™ S, SPECTROWORKS™ TOGETHER WITH UV-VIS SPECTROPHOTOMETER  
Date of issue: 30-Jul-2021  
Supersedes: -

Written by: Wenbo Wang

Wenbo Wang, QC

Date: 30 JULI 2021


Approved by: Søren Lund Kristensen

Søren Lund Kristensen, QA

Date: 30 JULI 2021

## POLYSTYRENE BEADS AS AN INDUSTRY STANDARD REFERENCE SYSTEM

Validated by external partner






## TEST AND VALIDATION IN REAL LIFE SETTINGS

- GLP Certified Particle analytical from Denmark has tested and validated the use of micro-NanoCuvette™ S and SpectroWorks™ using VWR UV6300-PC instrument.
- DLS was used for the reference instruments. DLS measurements were done by particle analytical using malvern zetasizer instrument.
- Mean particle size and concentration measurements of 6 different polystyrene latex beads were replicated for both of the instruments.



# PARTICLE SIZE DISTRIBUTION I – NANOCUVETTE™ S COMPARED TO DLS





| Size of Polystyrene latex beads | PDI DLS |       | Mean particle size (D50, nm) from DLS (% deviation from standard) |              | Mean particle size (nm) NanoCuvette™ S and SpectroWorks™ together with UV-vis spectrophotometer (% deviation from standard) |                |
|---------------------------------|---------|-------|---|--------------|---|----------------|
|                                 | Run 1   | Run 2 | Run 1   | Run 2        | Run 1   | Run 2          |
| 100 nm                          | 0.026   | 0.058 | 109 (9.0%)  | 107 (7.0%)   | 115.00 (15.0%)  | 148.90 (48.9%) |
| 460 nm                          | 0.025   | 0.048 | 442 (3.9%)  | 431 (6.3%)   | 414.10 (10.0%)  | 413.04 (10.2%) |
| 600 nm                          | 0.040   | 0.086 | 662 (10.3%)   | 633 (5.5%)   | 546.26 (9.0%)   | 549.85 (8.4%)  |
| 800 nm                          | 0.070   | 0.102 | 704 (12.0%)   | 730 (8.8%)   | 633.44 (20.8%)  | 637.50 (20.3%) |
| 1100 nm                         | 0.147   | 0.040 | 873 (20.6%)   | 953 (13.4%)  | 1097.66 (0.2%)  | 1083.85 (1.5%) |
| 3000 nm                         | 0.365   | 0.334 | 1690 (43.7%)  | 1480 (50.7%) | 3637.24 (21.4%)   | 3050.95 (1.9%) |

Table 1. Mean particle size (Vol.) of polystyrene latex beads measured by DLS and NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer


\*Data provided by Particle Analytical

## PARTICLE SIZE DISTRIBUTION II – NANOCUVETTE™ S COMPARED TO DLS



## PARTICLE CONCENTRATION ANALYSIS I – NANOCUVETTE™ S AND SPECTROWORKS™

Polystyrene latex beads concentration measurements





| Sample name                     | Concentration (Vol.)<br>Diluted from the stock | Concentration (Vol.)<br>From NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer. Result (% deviation) |                 |
|---------------------------------|--|--|-----------------|
|                                 |  | Run 1  | Run 2           |
| Polystyrene Latex beads 100 nm  | 0.01%  | 0.0099% (1.0%)   | 0.0081% (19.0%) |
| Polystyrene Latex beads 460 nm  | 0.003%   | 0.0038% (26.7%)  | 0.0038% (26.7%) |
| Polystyrene Latex beads 600 nm  | 0.003%   | 0.0026% (13.3%)  | 0.0026% (13.3%) |
| Polystyrene Latex beads 800 nm  | 0.003%   | 0.0030% (0.0%)   | 0.0031% (3.3%)  |
| Polystyrene Latex beads 1100nm  | 0.001%   | 0.0016% (60.0%)  | 0.0016% (60.0%) |
| Polystyrene Latex beads 3000 nm | 0.001%   | 0.0008% (20.0%)  | 0.0007% (30.0%) |

Table 2. Concentration of the diluted polystyrene latex beads obtained using NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer.

\*Data provided by Particle Analytical

## PARTICLE CONCENTRATION ANALYSIS II – NANOCUVETTE™ S AND SPECTROWORKS™



# PARTICLE SIZE ANALYSIS I - DYNAMIC LIGHT SCATTERING

## CERTIFICATE OF ANALYSIS

Customer: CphNano  
Material tested: Polystyrene Latex beads 100nm  
Batch: 2021/07/23  
Internal number: 32027

Analytical technique: Dynamic Light Scattering  
Method of analysis: Analyseplan  
Internal quality level: GMP

|               | D <sub>10%</sub> (nm)** | D <sub>50%</sub> (nm)** | D <sub>90%</sub> (nm)** | Z-average* | PDI*  |
|---------------|-------------------------|-------------------------|-------------------------|------------|-------|
| Run 1         | 79.4                    | 109                     | 159                     | 118        | 0.026 |
| Run 2         | 74.7                    | 107                     | 163                     | 118.1      | 0.058 |
| Average       | 76.9                    | 108                     | 161                     | 118.1      | 0.042 |
| Specification | -                       | -                       | -                       | -          | -     |
| Evaluation    | -                       | -                       | -                       | -          | -     |

\*The parameter is based on the intensity size distribution

\*\*The parameter is based on the volume size distribution

Written by: *Andi Guep* Date: 28 JULI 2021  
Reviewed by: 28 JULI 2021<sup>①</sup> *Wenbo Way* Date: 28 JULI 2021  
Approved by: 28 JULI 2021<sup>①</sup> *Wenbo Way* Date: 28 JULI 2021  
QC

## CERTIFICATE OF ANALYSIS

Customer: CphNano  
Material tested: Polystyrene Latex beads 460nm  
Batch: 2021/07/23  
Internal number: 32028

Analytical technique: Dynamic Light Scattering  
Method of analysis: Analyseplan  
Internal quality level: GMP

|               | D <sub>10%</sub> (nm)** | D <sub>50%</sub> (nm)** | D <sub>90%</sub> (nm)** | Z-average* | PDI*  |
|---------------|-------------------------|-------------------------|-------------------------|------------|-------|
| Run 1         | 320                     | 442                     | 605                     | 422.2      | 0.025 |
| Run 2         | 316                     | 431                     | 587                     | 413.6      | 0.048 |
| Average       | 318                     | 436                     | 597                     | 417.9      | 0.037 |
| Specification | -                       | -                       | -                       | -          | -     |
| Evaluation    | -                       | -                       | -                       | -          | -     |

\*The parameter is based on the intensity size distribution

\*\*The parameter is based on the volume size distribution

Written by: *Andi Guep* Date: 28 JULI 2021  
Reviewed by: *Wenbo Way* Date: 28 JULI 2021  
Approved by: *Wenbo Way* Date: 28 JULI 2021  
QC

**Conclusion: NanoCuvette™ S with developed software outperforms DLS.**

# PARTICLE SIZE ANALYSIS II - DYNAMIC LIGHT SCATTERING

## CERTIFICATE OF ANALYSIS

Customer: CphNano  
 Material tested: Polystyrene Latex beads 600nm  
 Batch: 2021/07/23  
 Internal number: 32029

Analytical technique: Dynamic Light Scattering  
 Method of analysis: Analyseplan  
 Internal quality level: GMP

|               | D <sub>10%</sub> (nm)** | D <sub>50%</sub> (nm)** | D <sub>90%</sub> (nm)** | Z-average* | PDI*  |
|---------------|-------------------------|-------------------------|-------------------------|------------|-------|
| Run 1         | 449                     | 662                     | 948                     | 569.5      | 0.04  |
| Run 2         | 470                     | 633                     | 831                     | 596.4      | 0.086 |
| Average       | 462                     | 646                     | 899                     | 582.9      | 0.063 |
| Specification | -                       | -                       | -                       | -          | -     |
| Evaluation    | -                       | -                       | -                       | -          | -     |

\*The parameter is based on the intensity size distribution

\*\*The parameter is based on the volume size distribution

Written by: *Paula Duggan*

Date: 28 JULI 2021

Reviewed by: *Wenbo Wang*

Date: 28 JULI 2021

Approved by: *Wenbo Wang*  
QC

Date: 28 JULI 2021

## CERTIFICATE OF ANALYSIS

Customer: CphNano  
 Material tested: Polystyrene Latex beads 800nm  
 Batch: 2021/07/23  
 Internal number: 32030

Analytical technique: Dynamic Light Scattering  
 Method of analysis: Analyseplan  
 Internal quality level: GMP

|               | D <sub>10%</sub> (nm)** | D <sub>50%</sub> (nm)** | D <sub>90%</sub> (nm)** | Z-average* | PDI*  |
|---------------|-------------------------|-------------------------|-------------------------|------------|-------|
| Run 1         | 496                     | 704                     | 948                     | 621.9      | 0.07  |
| Run 2         | 552                     | 730                     | 944                     | 684.9      | 0.102 |
| Average       | 529                     | 718                     | 946                     | 653.4      | 0.086 |
| Specification | -                       | -                       | -                       | -          | -     |
| Evaluation    | -                       | -                       | -                       | -          | -     |

\*The parameter is based on the intensity size distribution

\*\*The parameter is based on the volume size distribution

Written by: *Paula Duggan*

Date: 28 JULI 2021

Reviewed by: *Wenbo Wang*

Date: 28 JULI 2021

Approved by: *Wenbo Wang*  
QC

Date: 28 JULI 2021

**Conclusion: NanoCuvette™ S with developed software outperforms DLS.**

# PARTICLE SIZE ANALYSIS III - DYNAMIC LIGHT SCATTERING

## CERTIFICATE OF ANALYSIS

Customer: CphNano  
 Material tested: Polystyrene Latex beads 3000nm  
 Batch: 2021/07/23  
 Internal number: 32032

Analytical technique: Dynamic Light Scattering  
 Method of analysis: Analyseplan  
 Internal quality level: GMP

|               | D <sub>10%</sub> (nm)** | D <sub>50%</sub> (nm)** | D <sub>90%</sub> (nm)** | Z-average* | PDI*  |
|---------------|-------------------------|-------------------------|-------------------------|------------|-------|
| Run 1         | 1040                    | 1690                    | 4030                    | 2698       | 0.365 |
| Run 2         | 1150                    | 1480                    | 2290                    | 2756       | 0.334 |
| Average       | 1110                    | 1550                    | 3570                    | 2727       | 0.349 |
| Specification | -                       | -                       | -                       | -          | -     |
| Evaluation    | -                       | -                       | -                       | -          | -     |

\*The parameter is based on the intensity size distribution

\*\*The parameter is based on the volume size distribution

Written by: Antoine Dupont

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang  
 QC

Date: 28 JULI 2021

## CERTIFICATE OF ANALYSIS

Customer: CphNano  
 Material tested: Polystyrene Latex beads 1100nm  
 Batch: 2021/07/23  
 Internal number: 32031

Analytical technique: Dynamic Light Scattering  
 Method of analysis: Analyseplan  
 Internal quality level: GMP

|               | D <sub>10%</sub> (nm)** | D <sub>50%</sub> (nm)** | D <sub>90%</sub> (nm)** | Z-average* | PDI*  |
|---------------|-------------------------|-------------------------|-------------------------|------------|-------|
| Run 1         | 527                     | 873                     | 1250                    | 789.6      | 0.147 |
| Run 2         | 736                     | 953                     | 1230                    | 873.5      | 0.04  |
| Average       | 620                     | 922                     | 1240                    | 831.5      | 0.093 |
| Specification | -                       | -                       | -                       | -          | -     |
| Evaluation    | -                       | -                       | -                       | -          | -     |

\*The parameter is based on the intensity size distribution

\*\*The parameter is based on the volume size distribution

Written by: Antoine Dupont

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang  
 QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor  
 Quality agreement not in place

**Conclusion: NanoCuvette™ S with developed software outperforms DLS.**



# PARTICLE SIZE ANALYSIS I - NANOCUVETTE™ S AND SPECTROWORKS™

## Results


|                         |                 |  |
|-------------------------|-----------------|--|
| Refractive Index:       | 1.33612794 nD   | Material tested: Polystyrene Latex beads<br>100nm, 0.01%<br>Analytical technique: NanoCuvette™S,<br>SpectroWorks™ with Spectrophotometer |
| Reference fit quality:  | 39.4927 %       |  |
| Sample fit quality:     | 31.7062 %       |  |
| Mean Particle Diameter: | 115.00 nm       |  |
| Particle Concentration: | 0.0099 % (Vol.) |  |

## Results

|                         |                 |   |
|-------------------------|-----------------|---|
| Refractive Index:       | 1.36491180 nD   | Material tested: Polystyrene Latex beads<br>460nm, 0.003%<br>Analytical technique: NanoCuvette™S,<br>SpectroWorks™ with Spectrophotometer |
| Reference fit quality:  | 35.5789 %       |   |
| Sample fit quality:     | 18.6321 %       |   |
| Mean Particle Diameter: | 414.10 nm       |   |
| Particle Concentration: | 0.0038 % (Vol.) |   |


## Plots

Spectra Size distribution



## Plots

Spectra Size distribution




# PARTICLE SIZE ANALYSIS I - NANOCUVETTE™ S AND SPECTROWORKS™

## Results

|                         |                 |   |
|-------------------------|-----------------|---|
| Refractive Index:       | 1.35779775 nD   | Material tested: Polystyrene Latex beads<br>600nm, 0.003%<br>Analytical technique: NanoCuvette™S,<br>SpectroWorks™ with Spectrophotometer |
| Reference fit quality:  | 54.7016 %       |   |
| Sample fit quality:     | 28.0487 %       |   |
| Mean Particle Diameter: | 549.26 nm       |   |
| Particle Concentration: | 0.0026 % (Vol.) |   |

## Plots

Spectra Size distribution




## Results

|                         |                 |   |
|-------------------------|-----------------|---|
| Refractive Index:       | 1.34305982 nD   | Material tested: Polystyrene Latex beads<br>800nm, 0.003%<br>Analytical technique: NanoCuvette™S,<br>SpectroWorks™ with Spectrophotometer |
| Reference fit quality:  | 55.9293 %       |   |
| Sample fit quality:     | 24.9135 %       |   |
| Mean Particle Diameter: | 633.44 nm       |   |
| Particle Concentration: | 0.0030 % (Vol.) |   |

## Plots

Spectra Size distribution



Screenshots from developed software (SpectroWorks™)

# PARTICLE SIZE ANALYSIS I - NANOCUVETTE™ S AND SPECTROWORKS™


## Results

● Refractive Index: N/A  
 ● Reference fit quality: N/A  
 ● Sample fit quality: N/A  
 ● Mean Particle Diameter: 1097.66 nm  
 ● Particle Concentration: 0.0016 % (Vol.)

Material tested: Polystyrene Latex beads  
 110nm, 0.001%  
 Analytical technique: NanoCuvette™S,  
 SpectroWorks™ with Spectrophotometer

## Plots

Spectra Size distribution




## Results

● Refractive Index: 1.57637491 nD  
 ● Reference fit quality: 35.7368 %  
 ● Sample fit quality: 59.9385 %  
 ● Mean Particle Diameter: 3050.95 nm  
 ● Particle Concentration: 0.0007 % (Vol.)

Material tested: Polystyrene Latex beads  
 3000nm, 0.001%  
 Analytical technique: NanoCuvette™S,  
 SpectroWorks™ with Spectrophotometer

## Plots

Spectra Size distribution



## Poll Question 5

**Do you use any cloud-based software for your lab analysis? Or do you use any programming language like python and MATLAB to interpret/analyze your data?**

Upgrade your spectrophotometer with software

**SpectroWorks™**

# THE FUTURE OF UV-VIS LABWORK IS HERE.




Drag-and-drop easy  
data processing



Analysis in seconds



One-click export for  
reporting



Optimize and  
automate your UV-Vis  
instrument workflow.

**Online now.**

**SPECTROWORKS.COM**

## SpectroWorks™

Upgrade your spectrophotometer  
with NanoCuvette™ S

# Conclusion

## WHAT'S IN IT FOR YOU?

---

- Reliable accurate results.
- Works with existing UV-Vis spectrophotometers in the lab.
- No download of software with always updated online software platform.
- Cuvettes can be reused - check cuvette fit quality in SpectroWorks™.
- No extensive expertise needed for operation.
- Less time wasted.
- Low up-front cost.

*NB! Buy NanoCuvette™ S every day for 12 years before you have spent the same money as a DLS instrument.*





## BENEFITS FROM VWR COLLABORATION

---

- VWR/Avantor's UV-Vis spectrophotometer now supports particle size and concentration measurements.
- You can order cuvettes with our lab supplies.
- Standard lab purchasement contracts can be used.
- Data sheet, technical briefs, instruction briefs and knowledge base publicly available (<https://knowledge.cphnano.com/>).
- Shipping across Europe.
- Fits into the existing laboratory routines.



## WHY CHOOSE US?



### Supports Your Hardware

NanoCuvette™ One works with your existing UV-Vis spectrophotometer.



### Designed For Humans

Easy to use with no special training required.  
Online support included.




### Accurate and Instant Results

Results are calculated based on spectral analysis  
in less than 5 seconds.




### Unique Traceability

First in the industry, each cuvette is marked for  
traceability and access to online calibration data.



**SUPERIOR TECHNOLOGY WITH NO UP-FRONT COST  
OR TRAINING FOR YOUR LAB**



NanoCuvette™ S



**Sangita Khatri**

Research Scientist  
sakha@cphnano.com




**Christopher Lüscher, PhD**

Interim CEO  
chrlu@cphnano.com

Achieve easy, fast and reliable particle  
size analysis with NanoCuvette™ S

**Thank you!**

# Conclusions



You will receive an email with the link of the recording

You can listen to the recording of this presentation whenever you wish.

On the home page on VWR.com, select Events section, then enter the Webinars subsection.

You will receive a certificate of attendance 24 hours after the session.

If you need any technical information or for any other topic please write to [webinar@avantorsciences.com](mailto:webinar@avantorsciences.com)

# Thank you Questions?

[Webinar@avantorsciences.com](mailto:Webinar@avantorsciences.com)

Get your standard pack of  
10 NanoCuvette™ S by  
ordering.